

## Technical Note 997-103\_TN006\_090130

**OmniPro-Accept**

**Module: File import/export**

### OmniPro-Accept ASCII Format

#### Description

This document describes the ASCII file format for files generated with the OmniPro-Accept software. Following typographic conventions are used in the syntax column of tables below:

- Space character (ASCII character # 32)
- [tab] ASCII character # 9
- NbrOfCurves* Name of variable length character string

Each line will start with a character coding the type of information in that line:

- 'i' Information on number of measurements dumped; separators between individual measurements.
- '#' Comments line. This sign could also appear anywhere inside another line and marks the rest of the line as comments.
- '%' Label information line. The sign is immediately followed by a three letter code (for type of information) and then, depending on this code, the information in question. (The coding is specific for each kind of information.)
- '!' Operator comments.
- '=' A measurement point.

Every line of the output is terminated with [CR] [LF].

The output is divided into blocks. First is the header block: (Space character is denoted by an '·' character)

| Line | Syntax                                                               | Comment                                                                                                                 |
|------|----------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------|
| i    | :MSR·[tab] <i>NbrOfCurves</i> [tab]·#·No.·of·<br>measurement·in·file | <i>NbrOfCurves</i> determines the number of curves in the file and thus how many data blocks there will be in the dump. |
| ii   | :SYS·BDS·0·#·#·Beam·Data·Scanner·System                              |                                                                                                                         |

Then there is one block of curve data for each curve.

| Line | Syntax                                         | Comment                                      |
|------|------------------------------------------------|----------------------------------------------|
| 1    | #                                              |                                              |
| 2    | #·RFA300·ASCII·Measurement·Dump·(·BDS·format·) |                                              |
| 3    | #                                              |                                              |
| 4    | #·Measurement·number·[tab] <i>Nbr</i>          | <i>Nbr</i> = order number of the curve       |
| 5    | #                                              |                                              |
| 6    | %VNR·1.0                                       | Indicates the format version of ASCII dumps. |

|    |                                      |                                                                                                                                                                                                                                                         |
|----|--------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 7  | %MOD·[tab]Mode                       | Mode = identifies how the curve was measured.<br>Possible values are:<br>'FLM·' (Film)<br>'RAT·' (Ratio (RelativeDose))<br>'ABS·' (AbsoluteDose)<br>'INT·' (Integrated)<br>'UDF·' (Undefined/Isodose)<br>' '                                            |
| 8  | %TYP·[tab]Type                       | Type = identifies type of curve. Possible values are:<br>'SCN·' (Scan)<br>'ISO·' (Isodose)<br>'UDF·' (Undefined)                                                                                                                                        |
| 9  | %SCN·[tab]ScanType                   | ScanType = identifies type of scan. Possible values are:<br>'DPT·' (DepthDose)<br>'PRO·' (Profile)<br>'MTX·' (Matrix)<br>'DIA·' (Diagonal)<br>'UDF·' (Undefined/Isodose)                                                                                |
| 10 | %FLD·[tab]DetectorType               | Possible values of DetectorType are:<br>'ION·' (IonChamber)<br>'SEM·' (Semiconductor)<br>'UDF·' (Undefined)                                                                                                                                             |
| 11 | %DAT·[tab]DateOfCreation             | DateOfCreation = 'MM-DD-YYYY'                                                                                                                                                                                                                           |
| 12 | %TIM·[tab]TimeOfCreation             | TimeOfCreation = 'HH:MM:SS'                                                                                                                                                                                                                             |
| 13 | %FSZ·[tab]FieldWidth[tab]FieldHeight | FieldWidth and FieldHeight in mm                                                                                                                                                                                                                        |
| 14 | %BMT·[tab]RadType[tab]Energy         | Possible values of RadType are:<br>'COB·' (Cobalt)<br>'PHO·' (Photons)<br>'ELE·' (Electrons)<br>'UDF·' (Undefined)<br>Energy is a right-justified, seven-character string with energy value in MV or MeV always with one decimal.<br>Example: '...26.5' |
| 15 | %SSD·[tab]SSD                        | SSD in mm                                                                                                                                                                                                                                               |
| 16 | %BUP·[tab]BuildUp                    | BuildUp in 0.1 mm                                                                                                                                                                                                                                       |
| 17 | %BRD·[tab]BeamReferenceDist          | BeamReferenceDist in mm                                                                                                                                                                                                                                 |
| 18 | %FSH·[tab]Shape                      | Shape = the field shape. Possible values are:<br>'-1' (Undefined)<br>'0' (Circular)<br>'1' (Rectangular)<br>'2' (Irregular)                                                                                                                             |
| 19 | %ASC·[tab]AccessoryNbr               | Accessory number                                                                                                                                                                                                                                        |
| 20 | %WEG·[tab]WedgeNbr                   | Wedge number (angle)                                                                                                                                                                                                                                    |
| 21 | %GPO·[tab]GantryAngle                | GantryAngle in degrees                                                                                                                                                                                                                                  |
| 22 | %CPO·[tab]CollimatorAngle            | CollimatorAngle in degrees                                                                                                                                                                                                                              |

|    |                                                                                    |                                                                                                                                                                                                                                                                                                                                                                                                               |
|----|------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 23 | %MEA·[tab]MeasurementType                                                          | Possible values of <i>MeasurementType</i> are:<br>'-1' (Undefined)<br>'0' (Absolute dose)<br>'1' (Open depth)<br>'2' (Open profile)<br>'4' (Wedge)<br>'5' (Wedge depth)<br>'6' (Wedge profile)                                                                                                                                                                                                                |
| 24 | %PRD·[tab]ProfileDepth                                                             | <i>ProfileDepth</i> in 0.1 mm                                                                                                                                                                                                                                                                                                                                                                                 |
| 25 | %PTS·[tab]NbrOfPoints                                                              | <i>NbrOfPoints</i> = number of curve data points                                                                                                                                                                                                                                                                                                                                                              |
| 26 | %STS·[tab]StartX[tab]StartY[tab]StartZ·#·Start·Scan·values·in·mm·<br>(·X·,·Y·,·Z·) | <i>StartX</i> , <i>StartY</i> and <i>StartZ</i> are right-justified, seven-character strings with start values for each axis in origin-relative coordinates, in mm and always with one decimal. Example: '-100.0'.                                                                                                                                                                                            |
| 27 | %EDS·[tab]EndX[tab]EndY[tab]EndZ·#·End·Scan·values·in·mm·<br>(·X·,·Y·,·Z·)         | <i>EndX</i> , <i>EndY</i> and <i>EndZ</i> are right-justified, seven-character strings with end values for each axis in origin-relative coordinates, in mm and always with one decimal. Example: '-100.0'.                                                                                                                                                                                                    |
| 28 | !·CommentsLine1                                                                    | Operator comments, sixty characters long.                                                                                                                                                                                                                                                                                                                                                                     |
| 29 | !·CommentsLine2                                                                    | Operator comments, sixty characters long.                                                                                                                                                                                                                                                                                                                                                                     |
| 30 | #                                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                               |
| 31 | #[tab]·X·Y·Z·Dose                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                               |
| 32 | #                                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                               |
| 33 | ·[tab]XPos[tab]YPos[tab]ZPos[tab]Dose                                              | <i>XPos</i> , <i>YPos</i> and <i>ZPos</i> are right-justified, seven-character strings with origin-relative coordinates of the data point, in mm and always with one decimal. Example: '10.0'. <i>Dose</i> is also a right-justified, seven-character string but with the normalized dose value of the data point, in percent and always with one decimal. This field is repeated for every curve data point. |
| -  | :EOM·#·End·of·Measurement                                                          | -                                                                                                                                                                                                                                                                                                                                                                                                             |

The last block indicates the end of file.

| Line | Syntax             | Comment                    |
|------|--------------------|----------------------------|
| -    | :EOF·#·End·of·File | Indicates the end of file. |

## Example

This is an example of a depth dose measurement, made with an ion chamber, in a 100x100 mm 6 MV photon field. SSD is 1000 mm. The ASCII output was made using version 4.3 of the RFA-300 software.

```

:MSR·[tab]1[tab]#·No·of·measurement·in·file
:SYS·BDS·0·#·Beam·Data·Scanner·system
#
#·RFA300·ASCII·Measurement·Dump·(·BDS·format·)
#
#·Measurement·number·[tab]1
#
%VNR·1.0
%MOD·[tab]RAT
%TYP·[tab]SCN
%SCN·[tab]DPT
%FLD·[tab]ION

```

```
%DAT [tab] 02-03-1988
%TIM [tab] 14:15:25
%FSZ [tab] 100 [tab] 100
%BMT [tab] PHO [tab] .....6.0
%SSD [tab] 1000
%BUP [tab] 13
%BRD [tab] 0
%FSH [tab] 1
%ASC [tab] 0
%WEG [tab] 0
%GPO [tab] 0
%CPD [tab] 0
%MEA [tab] 1
%PRD [tab] 0
%PTS [tab] 25
%STS [tab] .....0 [tab] .....0 [tab] .....0 [tab] ..# Start Scan values in mm (X,Y,Z)
%EDS [tab] .....0 [tab] .....0 [tab] ..300.0 [tab] ..# End Scan values in mm (X,Y,Z)
! PDD data from Med. Phys. 7, 720 (1980) .....
! .....
#
# [tab] ..X.....Y.....Z.....Dose
#
= [tab] .....0 [tab] .....0 [tab] ...10.0 [tab] ...99.7
= [tab] .....0 [tab] .....0 [tab] ...20.0 [tab] ...100.0
= [tab] .....0 [tab] .....0 [tab] ...30.0 [tab] ...96.1
= [tab] .....0 [tab] .....0 [tab] ...40.0 [tab] ...91.8
= [tab] .....0 [tab] .....0 [tab] ...50.0 [tab] ...88.1
= [tab] .....0 [tab] .....0 [tab] ...60.0 [tab] ...83.9
= [tab] .....0 [tab] .....0 [tab] ...70.0 [tab] ...79.8
= [tab] .....0 [tab] .....0 [tab] ...80.0 [tab] ...75.7
= [tab] .....0 [tab] .....0 [tab] ...90.0 [tab] ...71.7
= [tab] .....0 [tab] .....0 [tab] ...100.0 [tab] ...67.8
= [tab] .....0 [tab] .....0 [tab] ...110.0 [tab] ...64.4
= [tab] .....0 [tab] .....0 [tab] ...120.0 [tab] ...60.5
= [tab] .....0 [tab] .....0 [tab] ...130.0 [tab] ...57.5
= [tab] .....0 [tab] .....0 [tab] ...140.0 [tab] ...54.2
= [tab] .....0 [tab] .....0 [tab] ...150.0 [tab] ...51.0
= [tab] .....0 [tab] .....0 [tab] ...160.0 [tab] ...48.1
= [tab] .....0 [tab] .....0 [tab] ...170.0 [tab] ...45.2
= [tab] .....0 [tab] .....0 [tab] ...180.0 [tab] ...42.5
= [tab] .....0 [tab] .....0 [tab] ...190.0 [tab] ...40.5
= [tab] .....0 [tab] .....0 [tab] ...200.0 [tab] ...38.3
= [tab] .....0 [tab] .....0 [tab] ...220.0 [tab] ...34.2
= [tab] .....0 [tab] .....0 [tab] ...240.0 [tab] ...30.5
= [tab] .....0 [tab] .....0 [tab] ...260.0 [tab] ...26.8
= [tab] .....0 [tab] .....0 [tab] ...280.0 [tab] ...24.1
= [tab] .....0 [tab] .....0 [tab] ...300.0 [tab] ...21.4
:EOM ..# End of Measurement
:EOF ..# End of File
```



In case of any further questions, please contact us under:

IBA Dosimetry GmbH  
Service Department  
Bahnhofstrasse 5  
DE-90592 Schwarzenbruck  
Germany

Service hotline:  
phone: +49 9128 607 38  
fax: + 49 9128 607 10  
mailto: [service@iba-group.com](mailto:service@iba-group.com)